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Plasma Separating Tubes (PST)

| Top Color | Additives | Principle | Uses |
|------------|-----------------------------------|--|--|
| Lavender | EDTA | -The strongest anti-coagulant - Ca^{+2} chelating agent - To preserve blood cells components | - Hematology - Blood bank (ABO) - HbA1C (Glycosylated Hb) |
| Light Blue | Sodium Citrate | Ca^{+2} chelating agent | - PT: Prothrombin Time - PTT: Partial Thromboplastin Time (in case of unexplained bleeding and liver disease) |
| Green | Sodium Heparin or Lithium Heparin | Heparin binds to Thrombin and inhibits the second step in the coagulation cascade (Prothrombin \rightarrow Thrombin) Heparin Fibrinogen $\not\rightarrow$ Fibrin | Enzymes Hormones Electrolytes (Na^+ , K^+ , Mg^+ , Cl^-) |

| Top Color | Additives | Principle | Uses |
|------------|--|---|--|
| Black | Sodium Citrate | Ca ²⁺ chelating agent | ESR (Erythrocyte Sedimentation Rate) to test how much inflammation in the patient, unexplained fever, Arthritis, Autoimmune Disorder |
| Gray | -Sodium Fluoride → -Potassium Oxalate → | Glycolysis inhibitor Anti-Coagulant | Glucose tests |
| Royal Blue | Heparin Na-EDTA | Anti-Coagulant Tube should not be contaminated with metals | Toxicology Trace Elements and metals |
| Yellow | ACD (Acid-Citrate Dextrose) | Anti-Coagulant | DNA Studies Paternity Test HLA Tissue Typing (Human Leukocyte Antigen) The body used this protein to differentiate the self-cells from non-self cells |

Serum Separating Tubes (SST)

| Top Tubes | Additives | Principle | Uses |
|-----------|--|---|--|
| Red | ----- Sometimes it has gel or silicon at the bottom of tube to reduce hemolysis | Enhancing the formation of blood clot | Serology -Antibodies -Hormones -Drugs Virology Chemistry Blood cross matching before blood transfusion |
| Gold | ----- It has gel at the bottom of the tube to separate serum from the blood | Serum separating from the blood through the gel in the tube | Serology Chemistry |

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1- اديتا "EDTA" acid acetic Ethylene diamine tetra

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BIBLIOGRAFIA

1. Kaplan L.A. Glucose. Kaplan A et al. Clin Chem The C.V. Mosby Co. St. Louis, Toronto, Princeton 1984; 1032-1036.
2. Trinder P. Ann Clin Biochem 1969; 6: 24-33.
3. Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACCC Press, 1995.
4. Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AACCC 2001.
5. Burts A et al. Tietz Textbook of Clinical Chemistry, 3rd ed AACCC 1999.
6. Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed AACCC 1995.

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Ref: MI41011 Cont. R: 6 x 30mL

CE GLUCOSE-LQ


Glucose
GOD-POD. Líquido

Quantitative determination of glucose IVD

Store at 2-8°C

PRINCIPLE OF THE METHOD

Glucose oxidase (GOD) catalyses the oxidation of glucose to gluconic acid. The formed hydrogen peroxide (H₂O₂), is detected by a chromogenic oxygen acceptor, phenol, 4-aminophenazone (4-AP) in the presence of peroxidase (POD):

The intensity of the color formed is proportional to the glucose concentration in the sample^{1,2}.

CLINICAL SIGNIFICANCE

Glucose is a major source of energy for most cells of the body; insulin facilitates glucose entry into the cells. Diabetes is a disease manifested by hyperglycemia; patients with diabetes demonstrate an inability to produce insulin^{1,2,3}. Clinical diagnosis should not be made on a single test result; it should integrate clinical and other laboratory data.

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SIGNS OF REAGENT DETERIORATION:

- Presence of particles and turbidity.
- Blank absorbance (A) at 505 nm ≥ 0.32.

ADDITIONAL EQUIPMENT

- MINDRAY BS-120 / BS-200E Autoanalyzer.
- General laboratory equipment.

SAMPLES

Serum or plasma, free of hemolysis⁴.

Serum should be removed from the clot as quickly as possible.

Stability of the sample: Glucose in serum or plasma is stable at 2-8°C for 3 days.

REFERENCE VALUES⁴

Serum or plasma:
60 – 110 mg/dL ≈ 3.33 – 6.10 mmol/L

These values are for orientation purpose; each laboratory should establish its own reference range.

MINDRAY BS-120 / BS-200E APPLICATION

PARAMETERS

| | | | |
|-----------------|---------------------|-----------------|---------------------|
| Test | GLU / GLU | R1 | 300 / 300 |
| IP | ** | R2 | * |
| Full Name | GLU / GLU | Sample volume | 1 / 1 |
| Standard IP | | R1 Blank | |
| Reac. Type | Endpoint / Endpoint | Mixed Igt Blank | |
| Pr. Wavelength | 510 / 505 | Linearity Range | 0 mg/dL - 500 mg/dL |
| Sec. Wavelength | | Linearity Limit | * |
| Direction | Increase / Increase | Substrate Limit | * |
| Reac. Time | 1_33 / 0_33 | Factor | * |
| Incuba. Time | | Prozone check | * |
| Units | mg/dL / mg/dL | q1 | q2 |
| Precision | Interco / Interco | q3 | q4 |
| | | PC | Abn |

CALIBRATION (Cal. - 3x weekly)

| | |
|-------------------------|-------------------------------------|
| Rule | One-point Linear / Two-point Linear |
| Sensitivity | 1 / 1 |
| Replicates | 2 / 2 |
| Interval (days) | 0 / 0 |
| Difference Limit | |
| SD | |
| Blank Response | |
| Error Limit | |
| Correlation Coefficient | |

Blank parameter must be performed in order to get good results in CALIB screen from main menu. The blank calibration is stable until 35 days. After this period the blank parameter must be performed again in order to validate the calibration.

QUALITY CONTROL

Control sera and calibrators are recommended to monitor the performance of assay procedures: SPINROL H Calibrador, SPINROL H Normal and Pathologic (Ref. 1002011, 1002120 and 1002210).

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BIBLIOGRAPHY

1. Kaplan L.A. Glucose. Kaplan A et al. Clin Chem The C.V. Mosby Co. St. Louis, Toronto, Princeton 1984; 1032-1036.
2. Trinder P. Ann Clin Biochem 1969; 6: 24-33.
3. Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACCC Press, 1995.
4. Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AACCC 2001.
5. Burts A et al. Tietz Textbook of Clinical Chemistry, 3rd ed AACCC 1999.
6. Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed AACCC 1995.

PACKAGING

Ref: MI41011 Cont. R: 6 x 30 mL

ABO Blood group

ABO blood group procedure

Take 7 tube for direct and indirect, Rh (D)_ **most immunogenic** , control .

 **Direct cell group in cell suspension** ——— to detect Ag


Label the test tube

Prepare 5% cell suspension
(2drops of RBC + 19 drops of saline .)

Add reagent antisera

Add 2 drops of reagent (Anti A, Anti B, Anti AB)


Add 1 drop of A₁Cell , B cell .

Add one drop of 5% patient red cell suspension

Centrifugation for 30 seconds,

Read and record agglutination reaction

| <i>Bloodgroup</i> | <i>Cell Ag</i> | <i>Anti A</i> | <i>Anti B</i> | <i>Anti AB</i> | |
|-------------------|----------------|---------------|---------------|----------------|--|
| A | A | + | - | + |  |
| B | B | - | + | + |  |
| AB | A and B | + | + | + |  |
| O | No Ag | - | - | - |  |

| <u>Cases</u> | <u>A₁ cell</u> | <u>B cell</u> | <u>AB in serum</u> | <u>Blood group</u> | |
|--------------|---------------------------|---------------|-----------------------------|--------------------|---|
| 1 | - | + | B | B |  |
| 2 | + | - | A | A |  |
| 3 | - | - | No Ab | AB |  |
| 4 | + | + | <u>A Ab and B Ab</u> | O |  |

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Store at 2-8°C

PRINCIPLE OF THE METHOD

Albumin in the presence of bromocresol green at a slightly acid pH, produces a colour change of the indicator from yellow-green to green-blue. The intensity of the color formed is proportional to the albumin concentration in the sample^{1,2,3,4}.

CLINICAL SIGNIFICANCE

One of the most important serum proteins produced in the liver is albumin.

This molecule has an extraordinarily wide range of functions, including nutrition, maintenance of oncotic pressure and transport of Ca²⁺, bilirubin, free fatty acid, drugs and steroids.

Variation in albumin levels indicate liver diseases, malnutrition, skin lesions such as dermatitis and burns or dehydration^{1,2,3,4}.

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The reagent is ready to use.

STORAGE AND STABILITY

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C, protected from light and contaminations prevented during their use.
Do not use reagents over the expiration date.

Signs of reagent deterioration:

- Presence of particles and turbidity.
- Blank absorbance (A) at 600 nm > 0,40.

ADDITIONAL EQUIPMENT

- SPIN 800 Autoanalyzer.
- General laboratory equipment.

SAMPLES

Serum or plasma, free of hemolysis: Stability 1 month at 2-8°C or 1 week at 15-25°C.

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QUALITY CONTROL

Control sera and calibrators are recommended to monitor the performance of assay procedures: SPINROL H Calibrator, SPINROL H Normal and Pathologic (Ref. 1002011, 1002120 and 1002210).

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PERFORMANCE CHARACTERISTICS

Measuring range: From detection limit of 0,0349 g/dL to linearity limit of 6 g/dL.

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BIBLIOGRAPHY

1. Gender S. Uric acid. Kaplan A et al. Clin Chem The C.V. Mosby Co. St. Louis, Toronto, Princeton 1984; 1268-1273 and 425.
2. Rodkey F L. Clin Chem 1965; 11: 478-487.
3. Webster D. Clin Chem. 1974; Acta 53: 109-115.
4. Doumas BT. Clin Chem. 1971; Acta 31: 87-96.
5. Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACCC Press, 1995.
6. Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AACCC 2001.
7. Burts A et al. Tietz Textbook of Clinical Chemistry, 3rd ed AACCC 1999.
8. Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed AACCC 1995.

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GPT (ALT)

NADH. Kinetico UV. FCC rec. Liquid.

Quantitative determination of alanine aminotransferase GPT (ALT) IVD

Store at 2-8°C

PRINCIPLE OF THE METHOD

Alanine aminotransferase (ALT) or Glutamate pyruvate transaminase (GPT) catalyses the reversible transfer of an amino group from alanine to α -ketoglutarate forming glutamate and pyruvate.

The pyruvate produced is reduced to lactate by lactate dehydrogenase (LDH) and NADH:

The rate of decrease in concentration of NADH, measured photometrically, is proportional to the catalytic concentration of ALT present in the sample¹.

CLINICAL SIGNIFICANCE

The ALT is a cellular enzyme, found in highest concentration in liver and kidney.

High levels are observed in hepatic disease like hepatitis, diseases of muscles and traumas, its better application is in the diagnosis of the diseases of the liver.

When they are used in conjunction with AST aid in the diagnosis of infarcts in the myocardium, since the value of the ALT stays within the normal limits in the presence of elevated levels of AST^{1,4,5}.

Clinical diagnosis should not be made on a single test result; it should integrate clinical and other laboratory data.

REAGENTS

| | | |
|-----|-----------------------------|-------------|
| R 1 | TRIS pH 7.8 | 100 mmol/L |
| | Lactate dehydrogenase (LDH) | 1200 U/L |
| | L-Alanine | 500 mmol/L |
| R 2 | NADH | 0.18 mmol/L |
| | α -Ketoglutarate | 15 mmol/L |

PREPARATION

DUAL MODE: Ready to use.

MONO MODE: Pour reagent 2 content over reagent 1. Mix thoroughly avoiding foam forming and it will be ready to use (W/R). Working reagent (W/R) stability: 21 days at 2-8°C or 72 hours at room temperature (15-25°C).

STORAGE AND STABILITY

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C, protected from light and contaminations prevented during their use.

Do not use reagents over the expiration date.

Signs of reagent deterioration:

- Presence of particles and turbidity.
- Blank absorbance (A) at 340 nm < 1.00.

ADDITIONAL EQUIPMENT

- Spectrophotometer or colorimeter measuring at 340 nm.
- Thermostatic bath at 25°C, 30°C ó 37°C (± 0.1°C).
- Matched cuvettes 1.0 cm light path.
- General laboratory equipment.

SAMPLES

Serum or plasma¹: Stability 7 days at 2-8°C.

INTERFERENCES

Anticoagulants currently in use like heparin, EDTA, oxalate and fluoride do not affect the results. Hemolysis interferes with the assay¹.

A list of drugs and other interfering substances with ALT determination has been reported^{1,3}.

APPLICATION SPINLAB 180

| | | | |
|-----------------|---------|------------------|--------|
| Name | GPT | Ref. male low | 0 U/L |
| Abstr. Name | GPT | Ref. male high | 40 U/L |
| Mode | Kinetic | Ref. female low | 0 U/L |
| Wavelength | 340 nm | Ref. female high | 32 U/L |
| Units | U/L | Ref. Pac. Low | * |
| Decimals | 0 U/L | Ref. Pac. High | * |
| Low Conc. | 5 U/L | Print value low | * |
| High Conc. | 500 U/L | Print value high | * |
| Calibrator name | | Control 1 | * |
| Prozone check | No | Control 2 | * |
| | | Control 3 | * |
| | | Control factor | 1.000 |
| | | Control other | 0.000 |

| | | | |
|------------------|--------------|------------------|--------------|
| DUAL MODE | | MONO MODE | |
| Sample blank | No | Sample blank | No |
| R1 bottle (mL) | 25 mL | R1 bottle (mL) | 25 mL |
| normal volume | 240 μ L | normal volume | 300 μ L |
| non volume | 240 μ L | non volume | 300 μ L |
| Sample | | Sample | |
| normal volume | 30.0 μ L | normal volume | 30.0 μ L |
| non volume | 15.0 μ L | non volume | 15.0 μ L |
| R2 bottle (mL) | | R2 bottle (mL) | |
| normal volume | 60.0 μ L | normal volume | 60.0 μ L |
| non volume | 60.0 μ L | non volume | 60.0 μ L |
| Predicción | No | | |
| Slope blank | No | | |
| Delay, min. sec. | 50, 180 sec. | Delay, min. sec. | 51, 175 sec. |
| Linearity limit | 10.0 % | Linearity limit | 10.0 % |
| Factor | -1746 | Factor | -1746 |

| | | | |
|-------------------|-----------|-------------------|-----------|
| Reagent blank | No | Reagent blank | No |
| Low Absorbance | <1.00 Abs | Low Absorbance | <1.00 Abs |
| High Absorbance | 3.000 Abs | High Absorbance | 3.000 Abs |
| R. Abs. L. Limit | <1.00 Abs | R. Abs. L. Limit | <1.00 Abs |
| R. Abs. H. Limit | 3.000 Abs | R. Abs. H. Limit | 3.000 Abs |
| R. Abs. Deviation | 3.000 Abs | R. Abs. Deviation | 3.000 Abs |

QUALITY CONTROL

Control sera are recommended to monitor the performance of assay procedures: SPINCONTROL H Normal and Pathologic (Ref. 1002120 and 1002210).

If control values are found outside the defined range, check the instrument, reagents and technique for problems.

Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

REFERENCE VALUES^{4,5}

| | | | |
|-------|--------------|--------|--------|
| | 25°C | 30°C | 37°C |
| Men | up to 22 U/L | 29 U/L | 40 U/L |
| Women | up to 18 U/L | 22 U/L | 32 U/L |

Normal newborns have been reported to show a reference range of up to double the adult, attributed to the neonate's hepatocytes. These values decline to adult levels by approximately 3 months of age. These values are for orientation purpose; each laboratory should establish its own reference range.

NOTES

SPINREACT has instruction sheets for several automatic analyzers. Instructions for many of them are available on request.

BIBLIOGRAPHY

- Murray R. *Alanine aminotransferase*. Kaplan A et al. Clin Chem The C.V. Mosby Co. St Louis. Toronto. Princeton 1984; 1088-1090.
- Young DS. *Effects of drugs on Clinical Lab. Tests*, 4th ed AACCC Press, 1995.
- Young DS. *Effects of disease on Clinical Lab. Tests*, 4th ed AACCC 2001.
- Burtis A et al. *Tietz Textbook of Clinical Chemistry*, 3rd ed AACCC 1999.
- Tietz N W et al. *Clinical Guide to Laboratory Tests*, 3rd ed AACCC 1995.

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- Nako H K. High-density lipoprotein (HDL) cholesterol. Kaplan A et al. Clin Chem The C.V. Mosby Co. St Louis. Toronto. Pinecote 1984; 1207-1213 and 437.
- Grove T H. Effect of reagent pH on Determination of HDL Cholesterol by precipitation with Sodium Phosphotungstate-magnesium Clin Chem 1979; 25:560.
- US National Cholesterol Education Program of the National Institutes of Health.
- Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed. AACCC Press, 1995.
- Young DS. Effects of disease on Clinical Lab. Tests, 4th ed. AACCC 2001.
- Burris A. et al. Tietz Textbook of Clinical Chemistry, 3rd ed. AACCC 1999.
- Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed. AACCC 1995.

PRESENTACION

Ref: 1001095 Cuv. R: 4 x 5 mL

SPINREACT,S.A.S, S.A.U. Ctra.Santa Coloma, 7-C-17176 SANT ESTEVE DE BAS (G) SPAIN
Tel. +34 972 69 08 00 Fax +34 972 69 00 00 e-mail: spinreact@spinreact.com

HDLC-P

HDL Coolesterol P Precipitating reagent

HDL cholesterol precipitating reagent
IVD

Store at 2-8°C

PRINCIPLE OF THE METHOD

The very low density (VLDL) and low density (LDL) lipoproteins from serum or plasma are precipitated by phosphotungstate in the presence of magnesium ions. After centrifugation the supernatant contains high density lipoproteins (HDL). The HDL cholesterol fraction is determined using the total cholesterol enzymatic reagent^{1,2}.

CLINICAL SIGNIFICANCE

HDL particles carry cholesterol from the cells back to the liver. HDL is known as "good cholesterol" because high levels are thought to lower the risk of heart disease.
A low HDL cholesterol levels, is considered a greater heart disease risk^{1,2}.
Clinical diagnosis should not be made on a single test result; it should integrate clinical and other laboratory data.

REAGENTS

| | | |
|-------------------------------|----------------------|-----------|
| R | Phosphotungstic acid | 14 mmol/L |
| Precipitating Reagent | Magnesium chloride | 2 mmol/L |
| Optional STD ^{1,2,3} | HDL Aq. Prim. Std. | 50 mg/dL |
| Optional reagent | Coolesterol CHOD-POD | |

PRECAUTIONS

R: H314-Causes severe skin burns and eye damage.
Follow the precautionary statements given in MSDS and label of the product.

PREPARATION

The reagent is ready to use.

STORAGE AND STABILITY

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C, protected from light and contaminants prevented during their use. Do not use reagents over the expiration date.

Signs of product deterioration:
- Presence of particles and turbidity.

ADDITIONAL EQUIPMENT

- Spectrophotometer or colorimeter measuring at 505 nm (500-550).
- Matched cuvettes 1,0 cm light path.
- General laboratory equipment.

SAMPLES

Serum or plasma¹: Free of hemolysis. Removed from the blood clot as soon as possible.
Stability: HDL Cholesterol is stable for 7 days at 2-8°C.

PROCEDURE

Precipitation^{1,2,3}
1. Pipette into a centrifuge tube:

| | |
|-------------|-----|
| R (µL) | 100 |
| Sample (mL) | 1,0 |

- Mix well; allow to stand for 10 min at room temperature.
- Centrifuge at 4000 r.p.m. for 20 min or 2 min at 12000 r.p.m.
- Collect the supernatant and proceed it as a sample in the total cholesterol determination.

CALCULATIONS

Follow the Instructions of the total cholesterol insert.

Calculated LDL-cholesterol (Friedewald)

LDLc = Total cholesterol + HDLc • (TG/5)

BS12-12 31/06/16

QUALITY CONTROL

Follow the Coolesterol reagent instructions of use.

REFERENCE VALUES¹

| HDL-cholesterol: | Men | Women |
|------------------|-------------|-------------|
| Lower risk | > 55 mg/dL | > 65 mg/dL |
| Standard risk | 35-55 mg/dL | 45-65 mg/dL |
| Increased risk | < 35 mg/dL | < 45 mg/dL |

LDL-cholesterol:

Suspected above: 150 mg/dL
Increased above: 190 mg/dL
These values are for orientation purpose; each laboratory should establish its own reference range.

PERFORMANCE CHARACTERISTICS

Measuring range: From detection limit of 1,57 mg/dL to linearity limit of 275 mg/dL.
If the results obtained were greater than linearity limit, dilute the sample 1/2 with NaCl 9 g/L and multiply the result by 2.

Precision:

| | Intra-assay (n=20) | Inter-assay (n=20) |
|--------------|--------------------|--------------------|
| Mean (mg/dL) | 33,9 75,8 | 34,8 75,4 |
| SD | 0,85 0,89 | 1,25 1,95 |
| CV (%) | 2,51 1,18 | 3,60 2,59 |

Sensitivity: 1 mg/dL = 0,0015 A.

Accuracy: Results obtained using SPINREACT reagents (y) did not show systematic differences when compared with other commercial reagents (x).

The results obtained using 50 samples were the following:

Correlation coefficient (r²): 0,99.

Regression equation: y= 0,9944x - 1,2346.

The results of the performance characteristics depend on the analyzer used.

INTERFERENCES

No interferences were observed with triglycerides up to 4 g/L¹.
A list of drugs and other interfering substances with HDL cholesterol determination has been reported by Young et. al^{1,2}.

NOTES

- The Precipitation procedure can be also performed with the half of reagent and sample volume.
- Calibration with the aqueous Standard may cause a systematic error in automatic procedures. In these cases, it is recommended to use a serum Calibrator.
- SPINREACT has instruction sheets for several automatic analyzers. Instructions for many of them are available on request.

BIBLIOGRAPHY

- Nako H K. High-density lipoprotein (HDL) cholesterol. Kaplan A et al. Clin Chem The C.V. Mosby Co. St Louis. Toronto. Pinecote 1984; 1207-1213 and 437.
- Grove T H. Effect of reagent pH on Determination of HDL Cholesterol by precipitation with Sodium Phosphotungstate-magnesium Clin Chem 25:560, 1979.
- US National Cholesterol Education Program of the National Institutes of Health.
- Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed. AACCC Press, 1995.
- Young DS. Effects of disease on Clinical Lab. Tests, 4th ed. AACCC 2001.
- Burris A. et al. Tietz Textbook of Clinical Chemistry, 3rd ed. AACCC 1999.
- Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed. AACCC 1995.

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BIBLIOGRAFIA

1. Kaplan L.A. Glucose. Kaplan A et al. Clin Chem The C.V. Mosby Co. St Louis. Toronto. Princeton 1984; 1032-1036.
2. Trinder P. Ann Clin Biochem 1969; 6: 24-33.
3. Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACCC Press, 1995.
4. Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AACCC 2001.
5. Burts A et al. Tietz Textbook of Clinical Chemistry, 3rd ed AACCC 1999.
6. Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed AACCC 1995.

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Ref: MI41011 Cont. R: 6 x 30mL


CE **GLUCOSE-LQ**
Glucose
 GOD-POD. Líquido

Quantitative determination of glucose
 IVD
 Store at 2-8°C

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The intensity of the color formed is proportional to the glucose concentration in the sample.^{1,2}

CLINICAL SIGNIFICANCE

Glucose is a major source of energy for most cells of the body; insulin facilitates glucose entry into the cells. Diabetes is a disease manifested by hyperglycemia; patients with diabetes demonstrate an inability to produce insulin.^{1,3,4} Clinical diagnosis should not be made on a single test result; it should integrate clinical and other laboratory data.

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SIGNS OF REAGENT DETERIORATION:

- Presence of particles and turbidity.
- Blank absorbance (A) at 505 nm > 0,32.

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CALIBRATION (Cal - Ir. Manual)

| | |
|-----------------|-------------------------------------|
| Rule | One-point Linear / Two-point Linear |
| Sensitivity | 1 / 1 |
| Replicates | 2 / 2 |
| Interval (days) | 0 / 0 |

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BIBLIOGRAPHY

1. Kaplan L.A. Glucose. Kaplan A et al. Clin Chem The C.V. Mosby Co. St Louis. Toronto. Princeton 1984; 1032-1036.
2. Trinder P. Ann Clin Biochem 1969; 6: 24-33.
3. Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACCC Press, 1995.
4. Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AACCC 2001.
5. Burts A et al. Tietz Textbook of Clinical Chemistry, 3rd ed AACCC 1999.
6. Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed AACCC 1995.

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SPINREACT UREA-37
Urea 37
o-Phthalaldehyde 37°C. Colorimetric

Quantitative determination of urea
IVD
Store at 2-8°C

PRINCIPLE OF THE METHOD
Urea in the sample reacts with o-phthalaldehyde in acid medium forming a coloured complex that can be measured by spectrophotometry:
Urea + o-phthalaldehyde + H⁺ → Isoindoline
The intensity of the color formed is proportional to the urea concentration in the sample¹.

CLINICAL SIGNIFICANCE
Urea is the final result of the metabolism of proteins; it is formed in the liver from their destruction.
It can appear the urea elevated in blood (uremia) in diets with excess of proteins, renal diseases, heart failure, gastrointestinal hemorrhage, dehydration or renal obstruction^{1,4}.
Clinical diagnosis should not be made on a single test result; it should integrate clinical and other laboratory data.

| Blank | Standard | Sample |
|----------------------------------|----------|--------|
| R 1 (mL) | 1,0 | 1,0 |
| Standard ^{1000 μM} (μL) | -- | 25 |
| Sample (μL) | -- | 25 |
| Final acid | -- | -- |
| R 2 (mL) | 1,0 | 1,0 |

3. Mix and incubate 15 min at 37°C.
4. Read the absorbance (A) against the Blank.
Calculations:
$$\frac{(A \text{ Sample} - A \text{ Blank})}{(A \text{ Standard} - A \text{ Blank})} \times 50 (\text{Std. conc.}) = \text{mg/dL urea in the sample}$$

10 mg/dL urea BUN divided by 0,466 = 21 mg/dL urea = 0,36 mmol/L urea¹.
Conversion factor: mg/dL x 0,1665 = mmol/L.

QUALITY CONTROL
Control sera are recommended to monitor the performance of assay procedures: SPINREACT H Normal and Pathologic (Ref. 1002120 and 1002210). If control values are found outside the defined range, check the instrument, reagents and calibrator for problems.
Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

REFERENCE VALUES¹
Serum: 15-45 mg/dL (2,49-7,49 mmol/L)
Urine: 20-35 g/24 h.
These values are for orientation purpose; each laboratory should establish its own reference range.

PERFORMANCE CHARACTERISTICS
Measuring range: From detection limit of 0,70 mg/dL to linearity limit of 200 mg/dL.
If the results obtained were greater than linearity limit, dilute the sample 1/2 with NaCl 9 g/L and multiply the result by 2.
Precision:

| Mean (mg/dL) | Intra-assay (n=20) | Inter-assay (n=20) |
|--------------|--------------------|--------------------|
| SD | 43,2 | 145 |
| CV (%) | 3,49 | 0,76 |
| | 1,52 | 1,91 |

Sensitivity: 1 mg/dL = 0,00459 A.
Accuracy: Results obtained using SPINREACT reagents (y) did not show systematic differences when compared with other commercial reagents (x).
The results obtained using 50 samples were the following:
Correlation coefficient (r): 0,9884
Regression equation: y=0,975x + 0,6
The results of the performance characteristics depend on the analyzer used.

INTERFERENCES
It is recommended to use heparin as anticoagulant. Do not use ammonium salts or fluorides¹.
A list of drugs and other interfering substances with urea determination has been reported by Young et al.^{2,3}.

NOTES
1. UREA CAL: Proceed carefully with this product because due its nature it can get contaminated easily.
2. Glassware and distilled water must be free of ammonia and ammonium salts.
3. Calibration with the aqueous standard may cause a systematic error in automatic procedures. In these cases, it is recommended to use a serum Calibrator.
4. Use clean disposable pipette tips for its dispensation.
5. SPINREACT has instruction sheets for several automatic analyzers. Instructions for many of them are available on request.

BIBLIOGRAPHY
1. Kaplan A, Urea. Kaplan A et al. Clin Chem The C.V. Mosby Co. St. Louis, Toronto, Princeton 1984; 1257-1260 and 437 and 416.
2. Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACCC Press, 1995.
3. Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AACCC 2001.
4. Burts A et al. Text: Textbook of Clinical Chemistry, 3rd ed AACCC 1999.
5. Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed AACCC 1995.

PACKAGING
Ref. 1001323 R1: 1 x 50 mL, R2: 1 x 50 mL, CAL: 1 x 2 mL.
Ref. 1001325 R1: 1 x 250 mL, R2: 1 x 250 mL, CAL: 1 x 5 mL.
Ref. 1001326 R1: 1 x 100 mL, R2: 1 x 100 mL, CAL: 1 x 2 mL.

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BIBLIOGRAFÍA
1. Kaplan A, Urea. Kaplan A et al. Clin Chem The C.V. Mosby Co. St. Louis, Toronto, Princeton 1984; 1257-1260 and 437 and 416.
2. Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACCC Press, 1995.
3. Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AACCC 2001.
4. Burts A et al. Text: Textbook of Clinical Chemistry, 3rd ed AACCC 1999.
5. Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed AACCC 1995.

PRESENTACIÓN
Ref. 1001323 R1: 1 x 50 mL, R2: 1 x 50 mL, CAL: 1 x 2 mL.
Ref. 1001325 R1: 1 x 250 mL, R2: 1 x 250 mL, CAL: 1 x 5 mL.
Ref. 1001326 R1: 1 x 100 mL, R2: 1 x 100 mL, CAL: 1 x 2 mL.

CE UREA-37

SPINREACT S.A.S.U Ctra Santa Colomba, 7 E-17176 SANT ESTEVE DE BAS (G) SPAIN
Tel. +34 972 69 08 00 Fax +34 972 69 02 96. e-mail: spreact@spreact.com

URIC ACID-LQ

Uric acid


Uricase-POD, Líquido

Quantitative determination of uric acid IVD

Store at 2-8°C

PRINCIPLE OF THE METHOD

Uric acid is oxidized by uricase to allantoin and hydrogen peroxide (H₂O₂), which under the influence of POD, 4-aminophenazone (4-AP) and 2,4-Dichlorophenol sulfonate (DCPS) forms a red quinonimine compound:

The intensity of the red color formed is proportional to the uric acid concentration in the sample¹².

CLINICAL SIGNIFICANCE

Uric acid and its salts are end products of the purine metabolism. With progressive renal insufficiency, there is retention in blood of urea, creatinine and uric acid. Elevated uric acid level may be indicative of renal insufficiency and is commonly associated with gout¹⁴. Clinical diagnosis should not be made on a single test result; it should integrate clinical and other laboratory data.

REAGENTS

| R 1 | Phosphate pH 7,4 | 50 mmol/L |
|---------|-------------------------------------|-----------|
| Buffer | 2,4-Dichlorophenol sulfonate (DCPS) | 4 mmol/L |
| | Uricase | 60 U/L |
| R 2 | Peroxidase (POD) | 660 U/L |
| Enzymes | Ascorbate oxidase | 200 U/L |
| | 4-Aminophenazone (4-AP) | 1 mmol/L |

PREPARATION

All the reagents are ready to use.

STORAGE AND STABILITY

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C, protected from light and contaminations prevented during their use. Do not use reagents over the expiration date.

Signs of reagent deterioration:

- Presence of particles and turbidity.
- Blank absorbance (A) at 520 nm $\geq 0,16$.

ADDITIONAL EQUIPMENT

- SPIN640 / SPIN640Plus Autoanalyzer.
- General laboratory equipment.

SAMPLES

- Serum or plasma¹: Stability 3-5 days at 2-8°C or 6 months at -20°C.
- Urine (24 h)¹: Stability 4 days at 15-25°C, pH > 8. Dilute sample 1/50 in distilled water. Mix. Multiply results by 50 (dilution factor).
If urine is cloudy, warm the specimen to 60°C for 10 min to dissolve precipitated urates and uric acid. Do not refrigerate.

REFERENCE VALUES¹

Serum or plasma:

Women 2,5 - 6,8 mg/dL \approx 149 - 405 μ mol/L
Men 3,6 - 7,7 mg/dL \approx 214 - 458 μ mol/L
Urine: 250 - 750 mg/24 h \approx 1,49 - 4,5 mmol/24 h

These values are for orientation purposes; each laboratory should establish its own reference range.

QUALITY CONTROL

Control sera and calibrators are recommended to monitor the performance of assay procedures: SPINCONTROL H Calibrator, SPINCONTROL H Normal and Pathologic (Ref. 1002011, 1002120 and 1002210).

If control values are found outside the defined range, check the instrument, reagents and technique for problems. Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

BARCODED REAGENTS LOAD MUST BE PRECEDED OF A SPINREACT "DATABASE" COPY INTO THE ANALYZER SOFTWARE. IT IS AVAILABLE UNDER REQUEST TO SPINREACT.

MOB8545- 030517

SPIN640 APPLICATION

| TEST INFORMATION | REAGENT VOLUME | |
|---------------------|----------------|----------------------|
| IP | ** | Vol. R1 150 |
| Test | URIC | Vol. R2 150 |
| Full Name | Uric Acid | Vol. R3 |
| Standard # | 1 | Vol. R4 |
| SAMPLE VOLUME | | |
| Vol. Sample Stand. | 7,5 | Decimal 0,01 Slope 1 |
| Vol. Sample Intra. | | Unit mg/dL Inter. 0 |
| Vol. Sample Dec. | | |
| REACTION PARAMETERS | | |
| Reac. Type | End Point | Direction Increase |
| Pr. Wave. | 525 | Reagent Blank 0-0 |
| Sec. Wave. | | React. Time 76-77 |

SPIN640Plus APPLICATION

| TEST PARAMETERS | URIC | Na | ** |
|-----------------|-------------|--------------------|-------------|
| Test | URIC | Na | URIC |
| Full Name | URICACID | Na | URIC |
| Reac. Type | End Point | Direction Increase | |
| Pr. Wave. | 525 | Sec. Wave. | |
| Unit | mg/dL | Decimal | 0,01 |
| Reagent Blank | 0-0 | React. Time | R1 - R2 |
| Vol. Sample | 7,5 μ L | R1 | 150 μ L |
| Increase | R2 | R2 | 150 μ L |
| Decrease | R1 | R1 | |
| Sample Blank | R1 | | |

The Calibration is stable until 36 days. After this period the Calibration must be performed again in order to obtain good results.

PERFORMANCE CHARACTERISTICS

Measuring range: From detection limit of 0,01647 mg/dL to linearity limit of 40 mg/dL.

If the results obtained were greater than linearity limit, dilute the sample 1/2 with NaCl 9 g/L and multiply the result by 2.

Precision:

| | Intra-assay (n=20) | | Inter-assay (n=20) | |
|-------------|--------------------|-------|--------------------|-------|
| Mean (mg/L) | 4,46 | 10,37 | 4,71 | 11,02 |
| SD | 0,02 | 0,05 | 0,06 | 0,15 |
| CV (%) | 0,46 | 0,44 | 1,20 | 1,37 |

Sensitivity: 1 mg/dL = 0,0323 (A).

Accuracy: Results obtained using SPINREACT reagents (y) did not show systematic differences when compared with other commercial reagents (x). The results obtained using 50 samples were the following:

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NOTES

1. Calibration with the aqueous standard may cause a systematic error in automatic procedures. In these cases, it is recommended to use a serum Calibrator.
2. Use clean disposable pipette tips for its dispensation.

BIBLIOGRAPHY

1. Schultz A. Uric acid. Kaplan A et al. Clin Chem The C.V. Mosby Co. St. Louis, Toronto, Princeton 1984; 1261-1266 and 418.
2. Fossati P et al. Clin Chem 1980;28:227-231.
3. Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACCPress, 1995.
4. Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AACCPress, 1995.
5. Burts A et al. Tietz Textbook of Clinical Chemistry, 3rd ed AACCPress, 1999.
6. Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed AACCPress, 1995.

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- 1- *Guyton and Hall textbook of medical physiology* (11th edition).
- 2- Kidd ,E. (2007). *Dental Caries: The Disease and Its Clinical Management*. Wiley-Blackwell. ISBN 978-1-4051-3889-5.
- 3- Edgar ,M.; Dawes ,C.; O'Mullane ,D. (2004). *Saliva and Oral Health*. British Dental Association. ISBN 0-904588-87-4.
- 4- **"Characterization of the edible bird's nest the "Caviar of the East""**. doi:10.1016/j.foodres.2005.02.008.
- 5- Maton ,Anthea (1993). *Human Biology and Health*. Prentice Hall. ISBN 0-13-981176-1.
- 6- Manuel Ramos-Casals; Haralampos M. Moutsopoulos; John H. Stone. *Sjogren's syndrome: Diagnosis and Therapeutics*. Springer, 2011.
- 7- Nicki R.; Walker ,Brian R.; Ralston ,Stuart H. ,(2010). *Davidson's principles and practice of medicine*. Edinburgh: Churchill Livingstone/Elsevier.
- 8- Paldino ,Michael; Mogilner ,Alon Y.; Tenner ,Michael S. (December 2003). "Intracranial hypotension syndrome: a comprehensive.
- 9- Saunders NR ,Habgood MD ,Dziegielewska KM (1999). "Barrier mechanisms in the brain, I. Adult brain". *Clin. Exp. Pharmacol. Physiol*.
- 10- Hajdu SI (2003). **"A note from history: discovery of the cerebrospinal fluid"**. *Annals of Clinical and Laboratory Science*. **33** (3): 334–6.